

Biochemical Diagnostics

Detectabuse™ AU/NZ MIX

Urine Drug Test Controls

for Hand Held Devices



Top Reasons to use Detectabuse™ Controls:

- Human urine-based matrix simulates patient samples
- Recommended by the majority of hand held device manufacturers
- 12 drug constituents within one control vial
- Each Detectabuse™ control contains: THC, Benzoyllecgonine (COC), PCP, Morphine (OPI), Methamphetamine, Amphetamine, Secobarbital (BAR), Oxazepam (BZO), Methadone, Methaqualone, Propoxyphene, Nortriptyline (TCA).
- Stable for up to 3 years frozen, or refrigerated 2-8°C (unopened)
- Ready to use in 5 ml vial size - minimum wastage
- Use Detectabuse™ controls to train staff in the accurate use of hand held devices

AUSTRALIAN-NEW ZEALAND MIX CUTOFF LEVELS

TARGET VALUES	-50%	Cutoff*	+50%
Delta-9-THC-COOH	25	50	75
Benzoyllecgonine	150	300	450
Phencyclidine (PCP)	12.5	25	37.5
Morphine	150	300	450
d-Methamphetamine	150	300	450
d-Amphetamine	150	300	450
Secobarbital	150	300	450
Oxazepam	100	200	300
Methadone	150	300	450
Methaqualone	150	300	450
Propoxyphene	150	300	450
Nortriptyline	500	1000	1500

The Detectabuse Advantage:

Detectabuse AU/NZ MIX controls are formulated to work with all known hand held single use drug-screening devices. We also offer a complete line of controls for immunoassay and GC/MS analysers.



Distributed in New Zealand by:
 Diagnostic Bioserve Ltd
 E:diagnostic@xtra.co.nz
 P: 07 5422325
 W:www.dbl.co.nz
 Contact: Raewyn Jones



Diagnostic Bioserve